

## QUICK REFERENCE: HAZARDS OF ULTRAVIOLET RADIATION

### UV radiation from laboratory equipment is more concentrated than natural occurring UV!!!

#### What is Ultraviolet Radiation??


Ultraviolet light (UV) has two levels of radiation:

- Non-ionizing: ranges from 40-400 nanometers and is the most common form of UV radiation being used in biomedical and microbiological research laboratories.
- Ionizing: ranges from 100-280nm and more concentrated than natural occurring UV, which poses a greater threat to personnel. Can be emitted from some types of lab equipment.

#### What is the Exposure Hazard of UV?

Exposure to UV light can injure both the eye and skin.


- Photokeratitis - inflammation of the cornea (outer protective coating of the eye). Can with very brief exposure or just a flash of UV light.
- Erythema - sunburn of the skin. Can occur within a few seconds of exposure to UV. Prolonged exposure can cause premature aging and cancer of the skin.



#### LIMITS of UV for sterilization:

- Germicidal lamp has **limited penetrating power**.
- The dynamic air stream in a BSC decreases the efficacy of the UV.
- **UV light does not penetrate soil, dust, or solid objects.**
- **Intensity of lamp diminishes over time**, decreasing germicidal activity.
- **Humidity levels** above 70% decrease the germicidal effect of UV.
- **Temperatures** below 77°F reduce the output of the germicidal wavelength.

**EHR does not recommend use of UV for decontamination.**



Instruments that contain UV lights – use with caution!



Transilluminator



Crosslinker

### Safety Practices and Precautions

**Minimize eye and skin exposure:**

- Always wear personal protective equipment (PPE) such as gloves, face shields, and lab coats
- Never occupy BSC while UV lamp is activated.
- Do not work in a room where a UV light is active.
- Use Transilluminators **ONLY** with protective shield in place.
- Do not use Crosslinkers if the door safety interlocking mechanism is not working.

Additional information and references can be found in the biosafety manual and on the EHR website.