

BSC MythBusters: Does Heat Really Affect Protection?



Environments For Science™

Kara Held, Ph.D.
Science Director

BSC MythBusters

- There are a lot of rules, “guidelines”, rumors, and myths for using a Biosafety Cabinet.
- Which are true?
- Which are not?
- WHY?



BAKER

Environments For Science™


Biosafety Cabinets (BSCs)

- A ventilated enclosure for work with biohazard agents assigned to biosafety levels 1 through 4.
- Provides 3 types of **CONTAINMENT**:
 - Personnel protection
 - Product protection
 - Environmental protection
- All BSCs contain: 1+ HEPA filters & motor/blower



Class II Type A2 BSCs

- Personnel, Product and, Environmental Protection from particulates and aerosols.
- Minimum 100 fpm intake air
- Partial recirculation



Normal BSC Operation




1. Strong front intake air
2. Side suction slots
3. Downward HEPA filtered air
4. Smoke split to front and back
5. Momentum Air Curtain

BAKER

Environments For Science™

MYTH: Don't use heat sources in a BSC



- Theoretical reasoning:
 - The gas needed to light the flame would be dangerous flowing around the hot BSC motor.
 - Heat changes airflow dynamics within the BSC that are critical for safety.

Volatile Gas release

- Using our previous research (see Volatile Chem Use in a BSC on BakerCo.com), we can calculate how much propane can be released into a BSC safely without explosion *given there is no spark or excessive heat*.

$$ER(\text{propane}) = \frac{Q_i * MW * LEL * 473}{403 * SG * S_F * 100} = \frac{267 \text{ cfm} * 44.1 \text{ g/mol} * 2.1\% * 473}{403 * 1.52 \text{ g/mL} * 10 * 100}$$
$$= 19 \text{ mL/min}$$

$$ER(\text{natural gas}) = \frac{Q_i * MW * LEL * 473}{403 * SG * S_F * 100} = \frac{267 \text{ cfm} * 18 \text{ g/mol} * 1\% * 473}{403 * 0.56 \text{ g/mL} * 10 * 100}$$
$$= 10 \text{ mL/min}$$

- Release rate of propane from a tank to a Bunsen burner
= 0.12 mL/min



BAKER

Environments For Science™

***If there is no spark, flame, cracked gas tubing, leaky valves, etc. Still not recommended!**

Volatile Gas spontaneous ignition

- Auto Ignition temperature of propane = 504°C
- Maximum motor/blower allowable temperature = 150°C

BUSTED*

BAKER

Environments For Science™

*If there is no spark, flame, cracked gas tubing, leaky valves, etc. Still not recommended!

Gas Explosion

- If you are fan of the show, you know they do the theoretical exercise, then the practical test. Here's where we try to blow up a BSC using propane and natural gas.
- We did NOT do this.
- But here are some examples of explosions:



BAKER

Environments For Science™

***Flammable gasses in a BSC are NOT recommended!**


Airflow disruption

- Visually
 - Smoke
- BSC performance/Containment
 - Particle Counting
 - NSF International Standard 49 Biological Testing Criteria



Sources of Heat

- Four Commonly used Heat Sources:
 - Standard Bunsen Burner
 - High Heat Bunsen Burner
 - Bacti-Cinerator
 - Spirit Lamp



Heat Sources in a BSC


- Initial Problems:
 - Lighting a flame in a windy environment
 - Keeping it lit for the duration of the experiment



Airflow Disruptions

Upward flow of air





High Heat Bunsen Burner




Normal Operation





Environments For Science™



High Heat Bunsen Burner



Normal Operation



BAKER

Environments For Science™

Spirit Lamp



Bacti-Cinerator




BAKER

Environments For Science™

Contamination Experimentation


- Tested on a 6 ft and 4 ft Class II Type A2 BSC (SterilGARD e3)
- Split work surface into zones



Note: All experiments were done with an 8" sash opening.

Particle Counting


- Can the BSC maintain ISO Class 5?



Spirit Lamp	PASS ✓
Bacti-Cinerator	PASS ✓
Bunsen Burner	FAIL X
High Heat BB	FAIL X

- Taller flames affected the Momentum Air Curtain and intake air

HEPA filters are not immune to heat




BAKER

Environments For Science™

Biological Testing

- Microbiological aerosol testing in accordance with NSF International Standard 49 testing
- Locations A, B, C, D for Product and Personnel; B, C for Cross.



NSF 49 Pass/Fail Criteria:

Personnel

- ≤ 10 CFUs Impingers
- ≤ 5 CFUs Slit samplers
- Control > 300 CFUs

Product


- ≤ 5 CFUs
- Control > 300 CFUs

Cross

- ≤ 2 CFUs 14+” away

Biological Testing

- Bunsen Burner
 - 6ft BSC:
 - Personnel = **4/4 PASS**
 - Product = **3/4 FAIL**
 - B,C,D
 - Cross = **2/2 FAIL**
 - B,C
 - 4ft BSC:
 - Personnel = **3/4 FAIL**
 - B,C,D
 - Product = **4/4 FAIL***
 - A, B, C*, D
 - Cross = **2/2 FAIL**
 - B,C




BAKER

Environments For Science™

***Major Failure**

Biological Testing

- High Heat Bunsen Burner
 - 6ft BSC:
 - Personnel = **1/4 FAIL**
 - **C**
 - Product = **3/4 FAIL**
 - **A, C, D**
 - Cross = **2/2 FAIL**
 - **B, C**
 - 4ft BSC:
 - Personnel = **2/4 FAIL***
 - **A, C***
 - Product = **4/4 FAIL**
 - **A, B, C, D**
 - Cross = **2/2 FAIL**
 - **B, C**




BAKER

Environments For Science™

***Major Failure**


Biological Testing

- Bacti-Cinerator
 - 6ft BSC:
 - Personnel = **1/4 FAIL**
 - **D**
 - Product = **1/4 FAIL**
 - **A**
 - Cross = **1/2 FAIL***
 - **C***
 - 4ft BSC:
 - Personnel = **2/4 FAIL**
 - **B,D**
 - Product = **4/4 FAIL**
 - **A, B, C, D**
 - Cross = **2/2 FAIL**
 - **B, C**



Biological Testing

- Spirit Lamp
 - 6ft BSC:
 - Personnel = **4/4 PASS**
 - Product = **1/4 FAIL***
 - **C***
 - Cross = **1/2 FAIL**
 - **C**
 - 4ft BSC:
 - Personnel = **1/4 FAIL**
 - **D**
 - Product = **2/4 FAIL**
 - **A, B**
 - Cross = **2/2 FAIL***
 - **B*, C***



BAKER


Environments For Science™

*Close to passing

*Major Failure

Biological Testing

- Cross Contamination Microbiological Testing also allows you to determine a “safe zone” away from the burner.



- *Spoiler:* There isn't one.

Observations of Note

- Heat built up quickly within the BSC ($>10^{\circ}\text{C}$)
 - The worst results were seen after the burners had been on for a while
- Results were worse in the 4ft BSC vs. 6ft BSC
 - Smaller work area = less heat dissipation = greater airflow disruptions
- Aerosol generation was affected by burners
 - Flow out of the nebulizer could be seen fluctuating

Conclusions

- Heat sources in a BSC cause problems
- Major airflow disruptions
- Personnel, Product and Cross Contamination Protection are compromised
- Gas in a BSC is NOT recommended!
- WHOLE work area is affected



-
- Use a Burner alone on bench **SAFE***✓
 - Work in a BSC without a flame **SAFE***✓
 - Burner inside a BSC **NOT ALWAYS SAFE X**

BAKER


Environments For Science™

*Given following proper procedures

Alternative Solutions

- Instead of using a flame inside a BSC:
 1. Segregate work to flame work on a bench, other work in BSC
 2. Eliminate use of a flame, try touch heat sources
 3. Disposable loops, sterilized toothpicks, etc.
 4. Reassessment of procedures

Don't be this guy →



BAKER

Environments For Science™

Upcoming BSC MythBusters

- 2+ people in a BSC?
- How much is “overloaded”?



BAKER

Environments For Science™



Environments For Science™

Help spread the word of BSC myths, and send me the ones you want busted!

kheld@bakerco.com