


Topic	Undesirable Example	Proper Example
<p>VAPORIZERS</p> <p>Calibrate - (1-2 years)</p> <p>Filled in fume hood or use a keyed refill tool.</p> <p>Vaporizers must be secured by frame or equivalent.</p> <p>Vaporizers are anesthetic specific</p>		 <p>The image shows a VetFlo Traditional Anesthesia System. It consists of a silver and purple vaporizer unit mounted on a black frame. The frame has a vertical scale on the left side. The vaporizer unit has a label that reads "VetFlo™ Traditional Anesthesia System" and "Kent Scientific CORPORATION". Below this, it says "For Veterinary Use Only" and "Use Only Isoflurane". The unit is secured by the frame, which is a proper example of how vaporizers should be handled.</p>